



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Universal protective gloves, Kevlar, up to 250°C, pair

Special Price

€80.00 was

€120.00

**plus VAT &
Shipping**

Product Images



Description

Waterproof, heat resistant, cut protection and good chemical resistance!

- heat resistant up to +250°C,
short-term up to +450°C
- not flammable
- the outer skin is made of highly resistant para-aramid synthetic fiber
- the middle layer is water and chemical resistant due to nitrile
- inner lining is made of soft cotton

Additional Information

| | |
|---|---|
| No. | 2-4207 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072002442 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | yellow |
| Material | Cotton Nitrile Polyphenylene terephthalamide (PPTA, Kevlar®) |
| sterile | No |
| Features gloves | waterproof with special mechanical resistance with special chemical resistance with special protection against thermal risks latex-free extra strong breathable angeroughened |
| TBST MAX | 450 °C |
| Gender (clothes, shoes, gloves) | Unisex |
| Glove serves | Personal security |
| Gloves provide protection against | Heat Cuts Stabs Chemicals (limited) |
| Protection category (PPE) | II Protective equipment for defense against hazards (not CAT. I and III) |
| Protective clothing for the protection of | Hands/ arms |
| Type/ Shape Glove | 5-Finger |
| resistant against | Water |
| with special protection against cold | No |

