



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® lab coat for ladies, lapel collar,
blended fabric, size 40**

**€47.50
plus VAT &
Shipping**

Product Images



Description

You are guaranteed to feel comfortable in it!

Comfortable to wear ladies coat in BP-Comfort quality made of 50% cotton, 50% polyester. The cotton part is on the inside, the PES part on the outside.

- 105 cm long
- approx. 215 g/m²
- concealed snap fastener placket
- boil-proof, chlorine-resistant, suitable for industrial laundering, easy-care
- suitable for tumble dryer
- 1 breast pocket, 2 side pockets
- Back slit, back belt
- figure-accentuating cut

Further sizes (32 to 58) on request!

Additional Information

No.	4-1387
Manufacturer (Brand)	neoLab
EAN	4058072022075
Transport temperature	Room temperature
Color	White
Material	Blended fabric
Length	105 cm
Features clothing	air-permeable lint-free figure-hugging breathable
Size Ladies Germany	40
Size Ladies International	M
Size International (Unisex)	M
Gender (clothes, shoes, gloves)	Ladies
Collar shape	Lapel collar
Material weight	215 g/m ²
Sleeve end	without cuff or elastic

