



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® cleanroom gown, PP white, universal
size 105 cm long, 10 pcs/pack**

**€73.00
plus VAT &
Shipping**

Product Images



Description

Cleanroom gown made of strong non-woven fabric, with the opening in the front, to close with press studs, with 1 sewn inside pocket. Universal size.

The outer end of the sleeves is formed by an elastic band bordered with an overlock seam.

Seam width: 4 mm

Number of stitches: 18 +/- 4 on 5 cm.

Material: Non-woven fabric made of pure polypropylene, antistatically treated

Elastic: Latex-free, synthetic elastomer

Sewing thread: 100% polyester 17 - 18 tex

Dermatological compatibility: No intolerance reactions

Classification : Category I: simple protective equipment

Length : 105 cm (one size fits all)

Storage: Store in a dry place at room temperature away from direct sunlight

Additional Information

No.	7-3070
Manufacturer (Brand)	neoLab
EAN	4058072038571
Transport temperature	Room temperature
Color	White
Material	Polypropylene (PP) Nonwoven
Length	105 cm
Features clothing	suitable for cleanroom air-permeable for 1xuse lint-free antistatic
Gender (clothes, shoes, gloves)	Ladies Men Unisex
Size Ladies International	One size fits all
Size International (Unisex)	One size fits all
Collar shape	Stand-up collar
Sleeve end	with elastic

