



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® cryogenic work apron, wearable on  
both sides, 107 cm long**

**€265.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

The protective shield against cryogenic splashes and vapors.  
Safety for the body!

- completely asbestos free
- 100% water repellent material
- wearable on both sides
- adjustable straps at waist and neck
- resistant to most chemical substances
- antistatic
- machine washable
- suitable for clean rooms

Width: 61 cm each

## Additional Information

---

No.	2-4316
Manufacturer (Brand)	neoLab
EAN	4058072003142
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polyamide 6.6 (PA 6.6. Nylon®) Synthetic leather
Wide	61 cm
Length	107 cm
Features clothing	waterproof with special chemical resistance suitable for work with liquid nitrogen antistatic

