




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® inspection gown, blue, size M

Special Price

€30.20 was

€37.80

**plus VAT &
Shipping**



Klettverschluss



Arm Gummi



Description

Disposable gown made of antistatic treated non-woven fabric, with elasticated arms and Velcro fastener at the front for comfortable dressing, in blue. Generous cut for best freedom of movement. The gown does not stick to the body and does not slip up when walking, due to antistatic treatment of the non-woven fabric.

Additional Information

No.	8-2008
Manufacturer (Brand)	neoLab
EAN	4058072050610
Transport temperature	Room temperature
Color	blue
Material	Nonwoven Polypropylene (PP)
Length	105 cm
Features clothing	air-permeable for 1xuse breathable antistatic
Gender (clothes, shoes, gloves)	Unisex Ladies Men
Size International (Unisex)	M
Smock protects against	Biological hazards and infectious material
Protective clothing for the protection of	Body
Sleeve end	with elastic

