



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



ABEBA ESD safety clogs white, size 39/40, pair

Special Price

€48.00 was

€60.40

**plus VAT &
Shipping**

Product Images



Description

Autoclavable safety clogs in white, suitable for laboratories in the pharmaceutical industry, food or operating rooms.

- CE, EN ISO 20347:2007
- A, E, SRA
- autoclavable, sterilizable, disinfectable special material
- foldable heel strap
- autoclavable at 121°C
- washable at 60°C
- comfortable fit
- anti-slip outsole
- ESD-compatible according to DIN EN 61340

Additional Information

No.	2-0003
Manufacturer (Brand)	ABEBA
EAN	4058072047146
autoclavable at 121°C	Yes
Color	White
Features clothing	water repellent with replaceable insole with antibacterial inner material ESD certified
Shoe size	39 40
Features sole	antistatic slip
Gender (clothes, shoes, gloves)	Unisex Ladies Men
Interior material	Elastoplastic (TPE)
Material sole	Elastoplastic (TPE)
Upper material	Elastoplastic (TPE)
Type shoes	open back with foldable heel strap
standardized according to DIN standard (DIN EN ISO)	DIN EN ISO 20347:2012
standardized according to DIN standard (DIN EN)	DIN EN 61340

