



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® respirator half mask medium, EN 140, without filter

€36.20
plus VAT &
Shipping

Product Images



Description

A wide range of economy through modular design!

These respirators from 3M are equipped with a lightweight dual filter system. A bayonet connection is used to attach the different filters.

Skin-friendly, anti-allergic. Optimal weight distribution.

Economical due to low maintenance.

Supplied without filters

The mask and its filters are tested according to EN 140, 141, 143, 371 and 372 standards and are "CE" marked.

Accessories: Various filters (Order No. 2-2783 to 2-2787)

Additional Information

No.	2-2781
Manufacturer (Brand)	neolab
EAN	4058072069087
Transport temperature	Room temperature
Color	gray
Material	Butyl rubber (IIR)
Protective masks and hoods Features	reusable can be used with different filters with exhalation valve gas/particle filter available as accessory CE-marked
Size (small, medium, large)	medium
Type respirator	Half-face respirator (filter available as accessory)
standardized according to DIN standard (DIN EN)	DIN EN 140:1998-12

