



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Mouth and nose cover, blended fabric, 60°C washable


Special Price

€24.90 was

€49.60

**plus VAT &
Shipping**

Product Images



Description

Breathable neoLabLine® mouth-nose covering (fabric mask) with 2 elastic bands to hang behind the ears

The material of the mask is made of a blended fabric (polyester 65% and cotton 35%) and is washable at 60°C. The fabric carries the oeko-Tex label.

We expressly point out that this is not a medical product. It cannot be relied upon "...to protect you or others from transmission of SARS-CoV-2..." (Source: Federal Institute for Drugs and Medical Devices).

The mask can be used several times.

However, the following cleaning instructions must be observed:

Fabric masks are best cleaned immediately after wearing - at least **60 degrees in the washing machine**. Switch off the ECO program, as otherwise sufficiently high temperatures may not be reached. Alternatively, place the mask in boiling water in a pot on the stove for at least five minutes. If small cracks or damage are visible after washing, do not reuse the mask.

Sterilizing a mask in the microwave is not recommended, as many masks contain metal parts and must not be placed in the microwave, as they will catch fire there.

Additional Information

No.	2-0291
Manufacturer (Brand)	neoLabLine
EAN	4058072122430
autoclavable at 121°C	No
Color	White
Material	Blended fabric polyester (PEs) cotton
Form	Strips
Protective masks and hoods Features	reusable breathable foldable washable at 60°C
Size (small, medium, large)	medium
Type respirator	Mouth and nose covering

