



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® hand dispenser for 2 rolls up to 2.54 inner Ø

€10.90
plus VAT &
Shipping

Product Images



Description

Handy dispenser made of resistant plastic, especially suitable for neoTape labeling tapes.

Suitable for 2 pieces of neoTape with 12.7 m length
or all adhesive tapes with maximum 2.54 cm inner diameter.
With 2 spacers and plastic clip for locking.

Additional Information

No.	2-6250
Manufacturer (Brand)	neoLab
EAN	4058072007737
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Weight	0.06 kg
Dimensions dispenser (LxWxH)	45 x 115 x 90 mm
for width tape/ roll	up to 25 mm
for DM tape inside	25.4 mm

