



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Multitape Dispenser for rolls up to 7.62 cm inner Ø

€120.00
plus VAT &
Shipping

Product Images



Description

Handy dispenser made of resistant plastic, especially suitable for neoTape labeling tapes.

Order No. 2-6180

Suitable for neoTape with 55 m length
or all adhesive tapes with maximum 7.6 cm inner diameter.

Lx W x H: 25 x 18.4 x 15.2 cm

Additional Information

| | |
|------------------------------|--|
| No. | 2-6180 |
| Manufacturer (Brand) | neoLabLine |
| EAN | 4058072007379 |
| Transport temperature | Room temperature |
| Material | Polymethylmetacrylate (PMMA) (acrylic glass) |
| Weight | 1.12 kg |
| for width tape/ roll | up to 25 mm |
| Dimensions dispenser (LxWxH) | 250 x 184 x 152 mm |
| for DM tape inside | 76 mm |

