




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Multitape Dispenser for rolls up to 2.54 cm inner Ø

€86.20
plus VAT &
Shipping

Product Images



Description

Handy dispenser made of resistant plastic, especially suitable for neoTape labeling tapes.

Order No. 2-6251

Suitable for up to 15 pieces of neoTape with a length of 12.7 m
or all adhesive tapes with a maximum inner diameter of 2.54 cm.

Lx W x H: 25 x 11.4 x 8.6 cm

Additional Information

No.	2-6251
Manufacturer (Brand)	neoLab
EAN	4058072007744
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Weight	0.42 kg
for width tape/ roll	up to 25 mm
Dimensions dispenser (LxWxH)	250 x 114 x 86 mm
for DM tape inside	25.4 mm

