



neoLab® Multiple dispenser

Product Images



Description

Dispenser for multiple adhesive tape rolls with 76 mm inner diameter and 89 mm width.
Separators keep the rolls apart.
Signing platform 102 x 152 mm with serrated cutting edge.
Heavy base plate with four non-slip rubber studs.
Made of resin coated steel.

Dimensions: 26 x 15 x 11 cm
Weight: 2.3 kg

Additional Information

No.	6-1130
Manufacturer (Brand)	neoLab
EAN	4058072024819
Material	Stahlpulverbeschichtet
Weight	2.3 kg
for width tape/ roll	bis 89 mm
Dimensions dispenser (LxWxH)	150 x 260 x 110 mm
for DM tape inside	76 mm

