



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® LabelManager 210D labeling device

€119.90
plus VAT &
Shipping

Product Images



Description

Lightweight model for the laboratory bench, can also be used as a mobile unit.

Easy navigation, Access keys to special characters.

Up to 2 lines per label. 6 font sizes.

Font styles: normal, bold, italic, shadow, outline, vertical, mirrored.

Underlining, 8 borders. 10 continuous copies.

For ribbon widths 6, 9, 12 mm.

6.8 mm/second.

Batteries: 6 x LR6, Mignon AA (Order No. 2-3024)

Accessories: Type ribbons 9 mm wide (Order No. 3-1190 ff.)

Type ribbons 12 mm wide (Order No. 3-1178 ff.)

Additional Information

No.	3-1225
Manufacturer (Brand)	neoLab
EAN	4058072012571
Transport temperature	Room temperature
Wide	465 mm
Height	64 mm
Labeling and cutting equipment features	mobile
Length	155 mm
Depth	155 mm
Weight	0.43 kg
Energy supply	Battery Mains
for width tape/ roll	up to 13 mm

