



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® adhesive tape 12.5 mm wide, white,
roll 36 m**

**€21.00
plus VAT &
Shipping**

Product Images



Description

Adhesive tape for glass and plastic utensils, racks, tubes and others.

Adheres to all clean surfaces including non-stick surfaces, and can be peeled off without leaving residue.

Markings made with ballpoint pen or ink will not smear.

Resistant to moisture and solvents, resistant to repeated heating.

Temperature resistant to -73°C.

Additional Information

No.	6-2480
Manufacturer (Brand)	neoLab
EAN	4058072025953
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Acrylic
Temp MAX	121 °C
Temp MIN	-73 °C
Wide	12.5 mm
Length	36 m
DM inside	25.4 mm
Labels and tapes properties	waterproof with special chemical resistance writable
suitable for	Labeling/ Marking Fixing/ Connecting

