



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brady nylon tape, 12.7 mm wide, black on white,  
4.8 m/roll**

**€63.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

Nylon fabric label tape for the BMP™21 portable label printer.

Suitable for all surfaces in the laboratory, for sample vials, tubes, microtiter plates and slides.

Permanent adhesive. Hand labeling possible.

Printing color: black on white

Temperature range: -196°C to +90°C

**Can also be used for autoclaving to a limited extent:** for 3 cycles of 1 hour each at 121°C

Suitable for applications in liquid nitrogen (only in conjunction with plastic tubes).

## Additional Information

---

No.	1-2205
Manufacturer (Brand)	Brady
EAN	4058072034290
autoclavable at 121°C	No
Color	White
Material	Polyamide 6.6 (PA 6.6. Nylon®)
Wide	12.7 mm
Length	4.8 m
TBST MAX	90 °C
TBST MIN	-196 °C
Color lettering	black
for model labeler	BMP™21 mobile

