



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brady nylon tape, 19.05 mm wide, black on white, 4.8 m/roll

**€75.20
plus VAT &
Shipping**

Product Images



Description

Nylon fabric label tape for the BMP™21 portable label printer.

Suitable for all surfaces in the laboratory, for sample vials, tubes, microtiter plates and slides.

Permanent adhesive. Hand labeling possible.

Printing color: black on white

Temperature range: -196°C to +90°C

Can also be used for autoclaving to a limited extent: for 3 cycles of 1 hour each at 121°C

Suitable for applications in liquid nitrogen (only in conjunction with plastic tubes).

Additional Information

No.	1-2206
Manufacturer (Brand)	Brady
EAN	4058072034306
autoclavable at 121°C	No
Color	White
Material	Polyamide 6.6 (PA 6.6. Nylon®)
Wide	19 mm
Length	4.8 m
TBST MAX	90 °C
TBST MIN	-196 °C
Color lettering	black
for model labeler	BMP™21 mobile

