



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brady label tape polyester, chemical resistant, width: 12.70 mm

€70.80
plus VAT &
Shipping

Product Images



Description

Chemical resistant polyester label tape for the portable Label Printer Order. No.: 1-2200 BMP21-LAB. Suitable for labeling slides, resistant to ethanol, DMSO and briefly to xylene. Printing color: black on white Temperature range: -80° to +120°C Application range: freezer and autoclave

Additional Information

No.	1-2215
Manufacturer (Brand)	Brady
EAN	4058072034399
autoclavable at 121°C	Yes
Color	White
Material	Polyester (PEs)
Wide	12.7 mm
Length	6.4 m
TBST MAX	120 °C
Labels and tapes properties	with special chemical resistance
TBST MIN	-40 °C
Color lettering	black
for model labeler	BMP™21 mobile

