



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brady label tape polyester, chemical resistant, width: 19.05 mm

**€81.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Chemical resistant polyester label tape for the portable Label Printer Order. No.: 1-2200 BMP21-LAB. Suitable for labeling slides, resistant to ethanol, DMSO and briefly to xylene. Printing color: black on white Temperature range: -80° to +120°C Application range: freezer and autoclave

## Additional Information

---

No.	1-2216
Manufacturer (Brand)	Brady
EAN	4058072034405
autoclavable at 121°C	Yes
Color	White
Material	Polyester (PEs)
Wide	19 mm
Length	6.4 m
TBST MAX	120 °C
Labels and tapes properties	with special chemical resistance
TBST MIN	-40 °C
Color lettering	black
for model labeler	BMP™21 mobile

