



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brady label tape polyester, transparent, width:
12.70 mm**

**€63.70
plus VAT &
Shipping**

Product Images



Description

Transparent polyester label tape for the portable Label Printer BMP21-LAB Order. No. 1-2200. Suitable for labeling vials, flasks and dishes.

Printing color: black on transparent

Temperature range: -80° to + 120°C

Application area: freezer and autoclave

Additional Information

No.	1-2212
Manufacturer (Brand)	Brady
EAN	4058072034368
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polyester (PEs)
Wide	12.7 mm
Length	6.4 m
TBST MAX	120 °C
TBST MIN	-80 °C
Color lettering	black
for model labeler	BMP™21 mobile

