



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Self-laminating tape 9.53 mm wide,  
black on white, 6.4 m/roll**

**€75.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

Polyester label tape for BMP™21 portable label printer.  
Suitable for flat, level surfaces, sample vials, tubes and slides.  
Permanent adhesive. Hand labeling possible. Extremely thin material.

Print color: black on white

Temperature range: -196°C to 130°C

Suitable for liquid nitrogen applications when label is applied as a wrap-around label and overlaps.

## Additional Information

---

No.	1-2201
Manufacturer (Brand)	neoLab
EAN	4058072034252
autoclavable at 121°C	Yes
Color	White
Material	Polyester (PEs)
Wide	9.5 mm
Length	6.4 m
TBST MAX	130 °C
TBST MIN	-196 °C
Color lettering	black
for model labeler	BMP™21 mobile

