



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Self-laminating adhesive tape, 12.7 mm wide, black on white, 6.4 m/roll

€86.10  
plus VAT &  
Shipping

### Product Images



## Description

---

Polyester label tape for BMP™21 portable label printer.  
Suitable for flat, level surfaces, sample vials, tubes and slides.  
Permanent adhesive. Hand labeling possible. Extremely thin material.

Print color: black on white

Temperature range: -196°C to 130°C

Suitable for liquid nitrogen applications when label is applied as a wrap-around label and overlaps.

## Additional Information

---

No.	1-2202
Manufacturer (Brand)	neoLab
EAN	4058072034269
autoclavable at 121°C	Yes
Color	White
Material	Polyester (PEs)
Wide	12.7 mm
Length	6.4 m
Color lettering	black
TBST MAX	130 °C
TBST MIN	-196 °C
for model labeler	BMP™21 mobile

