



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Self-laminating tape, 30.48 mm wide,
black on white, 4.2 m/roll**

**€106.80
plus VAT &
Shipping**

Product Images



Description

Self-laminating vinyl label tape for BMP™21 portable label printer and other suitable label printers. Suitable for all surfaces in the laboratory, for sample vials, tubes, microtiter plates and slides. Durable acrylic adhesive, self-laminating. Hand labeling possible.

Printing color: black on white

Temperature range: -80°C to +70°C

Suitable for applications in liquid nitrogen, if the label is applied overlapping as a wrap-around label.

Can also be used for labeling frozen PCR tubes.

Additional Information

No.	1-2209
Manufacturer (Brand)	neoLab
EAN	4058072034337
autoclavable at 121°C	No
Color	White
Material	Polyvinyl chloride (PVC. non-specific)
Wide	30.5 mm
Length	4.2 m
TBST MAX	70 °C
Labels and tapes properties	self-laminating
TBST MIN	-80 °C
Color lettering	black
for model labeler	BMP™21 mobile

