



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Youtility labels, 100 pcs./pack

€22.61
plus VAT &
Shipping

Product Images





Description

DURAN® YOUTILITY labels

The labeling area marked on the screw cap fits perfectly with the self-adhesive DURAN® YOUTILITY labels. The YOUTILITY labels are specially developed for and matched to the DURAN® YOUTILITY bottle system. A re-sealable, water-resistant box contains a roll of 200 labels each.

The labels are made of a tear-resistant, matte white polyester coated with a special adhesive. The writable label surface is versatile for writing on (copiers, printers, ballpoint pens, markers, etc.).

Chemical resistance and outstanding temperature property (-40 to +150° C).

Labels with Class 1 autoclaving/steam sterilization indicator (ANSI/AAMI/ISO 11140-1:2005).

Labels are non-toxic and lead-free and can be easily removed after use without leaving any residue.

Additional Information

No.	SD-2826
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3039
EAN	4058072162436
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	36 mm
Color	whitelblue
Material	Polyester (PEs)
Wide	36 mm
Form	rectangular
Length	36 mm
TBST MAX	150 °C
TBST MIN	-40 °C
Format/ Packing	Roll
Dimensions	36 x 70 mm

