



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Testo Mini Indicators, +60 to +82°C

€35.00 As low
as
€31.90
plus VAT &
Shipping

Product Images



Description

The temperature mini-indicators are self-adhesive, temperature-sensitive foils that react with color change to specific temperature overshoots. They are ideal for temperature monitoring of products and processes where a specific temperature should not be exceeded, e.g. for moving or small objects or for long-term monitoring. The mini-indicators in use The temperature mini-indicators are supplied in a booklet of 10 each. They can be easily detached from the booklet like stickers and stuck onto the object to be measured. On the mini-indicators you can see the temperature points +60 °C; +65 °C; +71 °C; +77 °C; +82 °C. As soon as a certain temperature point is exceeded, the mini-indicator turns color in the relevant points within 2 to 3 seconds. The coloration is irreversible: once the temperature has been exceeded, the mini-indicator does not change color even if the temperature drops again. In this way, critical temperature overshoots can be read even after a longer period of time. Do you need mini-indicators for monitoring other temperature ranges? We also have the following measuring ranges in our product range: +88 ... +110 °C +116 ... +138 °C +143 ... +166 °C +171 ... +193 °C +199 ... +224 °C

Additional Information

No.	TO-0279
Manufacturer (Brand)	Testo
EAN	4066292351097
Manufacturer EAN	4029547006890
Gross weight	0.01
Transport temperature	Room temperature
Net weight	0.01000 kg
Country of origin	Great Britain
Packing width	0.10500 m
Packing height	0.00200 m
Packing depth	0.07000 m
Customs tariff number	38249993900

