



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® fabric adhesive tape, yellow, 19 mm
wide, 50 m/roll**

**€34.20
plus VAT &
Shipping**

Product Images



Description

A PVC-coated fabric adhesive tape for everyday use where objects are to be firmly and permanently joined, closed and marked.

Universally applicable.

Tear-resistant, yet very conformable.

Temperature resistant up to 90°C

Adhesive strength: 490 g/cm

Adhesive: Synthetic rubber

Accessories: Dispenser (Order No. 1-7501)

Additional Information

No.	2-3255
Manufacturer (Brand)	neoLab
EAN	4058072001124
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow
Material	Polyvinyl chloride (PVC. non-specific)
Temp MAX	90 °C
Wide	19 mm
Length	50 m
Labels and tapes properties	inscribable
Thickness	0.3 mm
Adhesive strength	0.26
suitable for	Fixing/ Connecting Labeling/ Marking

