



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® fabric tape tear-resistant, 50 mm, 50 m/roll

€25.60
plus VAT &
Shipping

Product Images



Description

Fabric tape tear-resistant

Suitable for sealing and covering large areas of metals and plastics, adheres securely and reliably. The tape can be easily labeled with a pencil or ballpoint pen .

Roll: 50 m long, 50 mm wide

Accessories: dispenser for large tapes (Order No. 2-6266)

Additional Information

No.	2-6263
Manufacturer (Brand)	neoLab
EAN	4058072007799
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	light gray
Material	Polyester (PEs) acrylate
Wide	50 mm
Length	50 m
DM inside	75 mm
Labels and tapes properties	inscribable
suitable for	Fixing/ Connecting Labeling/ Marking

