




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bohlender PTFE sealing tape, 12 mm wide, 0.1 mm thick, 12 m/roll

**€2.45** As low as  
**€1.56**  
plus VAT &  
Shipping

### Product Images



## Description

---

Made of pure PTFE for all threads made of steel, copper, brass, plastic, etc.

The material is resistant to almost all liquids, even the most aggressive media, does not swell, stick and does not age. Screw connections can thus be loosened without problems even after years.

Temperature resistance:

-190°C to +250°C

## Additional Information

---

No.	ER-1721
Manufacturer (Brand)	Bohlender
old neoLab article no.	2-6259
EAN	4058072246068
Transport temperature	Room temperature
Color	White
Material	Polytetrafluoroethylene (PTFE)
Wide	12 mm
Length	12 mm
suitable for	variable applications
TBST MAX	250 °C
TBST MIN	-190 °C
Labels and tapes properties	with special chemical resistance
Thickness	0.1 mm

