



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Temperature-resistant adhesive tape 25 mm, white, 7.6 m/roll

€50.90
plus VAT &
Shipping

Product Images



Description

High quality white textile based adhesive tape with Special adhesive, self-adhesive . Ideal product for extreme temperatures. Suitable for cryogenic applications as well as for dry incubators, sterilizers, etc. Easy to mark with pencil, ballpoint pen, felt-tip pen, etc.; various widths available.

Additional Information

| | |
|-----------------------------|--|
| No. | 2-6262 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072007782 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | White |
| Material | Cellulose mixed ester glass fiber reinforced |
| Temp MAX | 287 °C |
| Temp MIN | -231 °C |
| Wide | 25 mm |
| Length | 7.6 m |
| Labels and tapes properties | inscribable |
| suitable for | Labeling/ Marking |

