



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PARAFILM M, 10 cm wide, 38 m per roll

€42.50
plus VAT &
Shipping

Product Images



Description

PARAFILM® M sealing film can be stretched up to 200 % and also hugs irregular shapes and surfaces.

Without plasticizers.

Resistant to salt solutions, inorganic acids and alkalis for up to 48 hours.

Continuous service temperature -45°C to +50°C.

Additional Information

No.	3-1011
Manufacturer (Brand)	neoLab
EAN	4058072011901
Transport temperature	Room temperature
Color	Blue White
Material	Paraffin Low Density Polyethylene (LDPE)
Temp MAX	50 °C
Temp MIN	-45 °C
Wide	10 cm
Length	38 m
Labels and tapes properties	self-adhesive with special chemical resistance stretchable
Format/ Packing	Roll
Type foil	Closure material

