



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® film sealing device incl. 1 roll, 29 cm
sealing width**

**€93.80
plus VAT &
Shipping**

Product Images



Description

Welds and separates in a single operation.

With powerful air extraction, therefore also ideal for packaging all types of frozen goods.

Control light for sealing process with automatic switch-off.

Unit is supplied complete with a roll of film.

Sealing seam width: 29 cm

Connection: 230 Volt/60 Hz

Additional Information

No.	2-1095
Manufacturer (Brand)	neoLab
EAN	4058072055486
Color	White
Wide	100 mm
Length	290 mm
Weight	2 kg
Power	160 W
Refill package/article	2-1096
Welding seam length	290 mm
Voltage	230 V
for material foil	Low-density polyethylene (LDPE)

