




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Felt pen with varnish color, green

€12.10
plus VAT &
Shipping

Product Images



Description

Felt marker with fast drying lacquer paint, for marking metal, wood, plastic. Glass, rubber, cardboard. Adheres to rough, smooth, damp, oily and dry surfaces. The ink does not fade, does not rub off and is heat and weather resistant. The felt tip is replaceable.

Additional Information

No.	2-5071
Manufacturer (Brand)	neoLab
EAN	4058072004187
Transport temperature	Room temperature
Color	green
Marker properties	for wet surfaces for glass
Spare part	2-5073
Line width	3.2 mm
Type pin	Felt-tip pen

