



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® 3000 permanent marker, set of 10, line width 1.5-3mm

€29.90
plus VAT &
Shipping

Product Images



Description

The **edding permanent markers** with long-lasting and abrasion-resistant ink prove to be a reliable all-in-one solution for permanent marking and labeling in industry, laboratories and offices.

The waterproof and low-odor ink is suitable for all common metal, plastic and glass surfaces

- fast drying
- various colors and line widths

Additional Information

No.	1-1981
Manufacturer (Brand)	neoLab
EAN	4058072156664
Color	yellow orange red pink violet light blue blue green brown black
Marker properties	in set for glass
Number of parts in set	10
Line width	3 mm
Type pin	Felt-tip pen

