




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Blotting Membrane CN, 0.22 µm, 3 m/roll

€405.00
plus VAT &
Shipping

Product Images



Description

Blotting membranes are ideal for a wide range of blotting methods such as electrotransfer, semidry and capillary blotting.

Easy to cut to the desired gel dimension.

Optimal for Western, Southern and Northern blotting.

High binding capacity, low background noise.

Additional Information

No.	E-0026
Manufacturer (Brand)	neoLab
EAN	4058072062057
Transport temperature	Room temperature
Wide	30 cm
Length	300 cm
Material	Cellulose nitrate
Material membrane	Cellulose nitrate
Pore size	0.22 µm
Format/ Packing	Roll
Type filter	Blotting membrane Transfer membrane

