



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**qpore® transfer membrane nitrocellulose,
100x100 mm, 0.22 µm, 100 pcs/pack**

**€319.00
plus VAT &
Shipping**

Product Images



Description

Transfer membranes made of CN are suitable for a wide range of applications: Western, Northern and Southern blotting, protein dot blots and nucleic acid detection. The transfer membranes are very robust and hydrophobic. They increase sensitivity due to low background. Practical in 10x10 cm format

Additional Information

No.	6-0049
Manufacturer (Brand)	qpore
EAN	4058072169060
Transport temperature	Room temperature
Wide	10 cm
Length	10 cm
Material	Nitrocellulose
Format/ Packing	Sheet/ Sheet
Material membrane	Nitrocellulose
Pore size	0.22 µm
Type filter	Transfer membrane

