



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Falcon® bottle top filter 500 ml, orange, CA 0.22
µm, HW 33, 12 pcs./pack**

**€138.31
plus VAT &
Shipping**

Product Images



Description

Bottle-top filters 500 ml can be used in conjunction with glass or polystyrene bottles. Two different designs, suitable for bottles with neck widths of 33 mm or 45 mm. Choice of different membrane types and pore sizes to suit the application.

For easy product identification, the integrated membrane type and pore size is printed on each system.

- all filters are sterile, certified and pyrogen-free
- Polystyrene refill container and bottle
- Filter diameter 70 mm
- Filter area: square

Membrane material:

CA - cellulose acetate

Additional Information

No.	XC-0715
Manufacturer (Brand)	Falcon
old neoLab article no.	7-4195
EAN	4058072444204
Transport temperature	Room temperature
Color	orange
Volume	500 ml
sterile	Yes
Material housing	Polystyrene (PS)
Material membrane	Cellulose acetate (CA)
Pore size	0.22 µm
Type connection	Hose olive
Type filter	Bottle top filter
for DM hose inside 1	8 - 12 mm

