



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Bottle-Top Filter System, 0.22 µm, 5-10 l, 10 pcs./pack

€247.00
plus VAT &
Shipping

Product Images



Description

Stericap™ PLUS for vacuum filtration of tissue culture media, microbiological media and biological solutions.
With Express® PLUS membrane (PES) for high filtration speed and high total throughput (5-10 liters)
For bottles with an inner neck diameter of 20 to 76 mm.
Individually sterile packed.

Additional Information

No.	MP-1019
Manufacturer (Brand)	neoLab
old neoLab article no.	7-9059
EAN	4058072335472
Transport temperature	Room temperature
Color	blue
Volume	10000 ml
sterile	Yes
Material housing	Polystyrene (PS)
Material membrane	Polyethersulfone (PES)
Pore size	0.22 µm
for DM hose inside 1	8 - 12 mm
Type connection	Hose olive
Type filter	Bottle top filter

