



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Millipore® Bottle-Top-Filter for 45 mm bottle neck, 250 ml sterile

€126.00
plus VAT &
Shipping

Product Images



Description

Suitable for filtration of tissue culture media, microbiological media and biological solutions.
High filtration speed at high total throughput with theExpress®PLUS membrane (PES).

Additional Information

| | |
|-----------------------|------------------------|
| No. | MP-1744 |
| Manufacturer (Brand) | Millipore |
| EAN | 4058072049515 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Color | blue |
| Volume | 250 ml |
| sterile | Yes |
| Material housing | Polystyrene (PS) |
| Material membrane | Polyethersulfone (PES) |
| Pore size | 0.22 µm |
| for DM hose inside 1 | 8 - 12 mm |
| Type connection | Hose olive |
| Type filter | Bottle top filter |
| for thread/DM | GL 45 |

