



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Bottle-Top-Filter for 45 mm bottle neck, 150 ml sterile, 12 pcs./pack

€129.00
plus VAT &
Shipping

Product Images



Description

Suitable for filtration of tissue culture media, microbiological media and biological solutions.
High filtration speed at high total throughput with theExpress@PLUS membrane (PES).

Additional Information

No.	MP-1023
Manufacturer (Brand)	neoLab
old neoLab article no.	7-9067
EAN	4058072335489
Transport temperature	Room temperature
Color	blue
Volume	150 ml
sterile	Yes
Material housing	Polystyrene (PS)
Material membrane	Polyethersulfone (PES)
Pore size	0.22 µm
Type connection	Hose olive
Type filter	Bottle top filter
for DM hose inside 1	8 - 12 mm
for thread/DM	GL 45

