



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Nalgene Rapid-Flow™ Sterile Disposable Bottle Top Filter with SFCA Membrane 45 mm

**€131.04**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES Membrane Perform sterile filtration of cell culture media and sera using Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES Membrane. Disposable Product: Yes Filter Type: Membrane Material (Membrane): SFCA Pore size: 0.2 µm Fits bottles with neck size (metric): 45 mm Number per case: 12 Number per package: 12 Type: Bottle top filter Volume (metric): 500 mL Diameter (metric) Membrane: 75 mm Unit Size: Case of 12

## Additional Information

---

No.	CD-0717
Manufacturer (Brand)	Nalgene
VGKL number	108729111
EAN	4058072219116
Gross weight	1.5
Transport temperature	Room temperature
Country of origin	Mexico
sterile	Yes
Packaging volume	0.02500000 ccm
Customs tariff number	84212980

