



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**qpore® Bottle-Top-Filter made of CA, sterile,
0.45 µm, 250 ml, 24 pcs/pack**

**€169.00
plus VAT &
Shipping**

Product Images



Description

qpore® offers a comprehensive range of bottle-top filters made of different materials, which are particularly suitable for the filtration of larger sample volumes. All membranes are manufactured under the highest quality standards from the best raw materials.

These **sterile qpore® bottle-top filters with cellulose acetate membrane (CA)** have fast flow rates and very low protein binding, making them ideally suited for filtration of cell culture media. With a GL 45 screw cap and tubing connection, these bottle-top filters are suitable for negative pressure filtration applications directly into the sterile bottle.

Features:

- Blue color coding
- Quick connector for tubing with inner Ø 8-12 mm
- Graduated filter top
- Sterile, pyrogen free, detergent free
- Bottle-top filters are individually sterile packaged in 24 pieces
- Packaging with cut open aid; thus can be opened without further aids

Fields of application:

- Sterile filtration of cell culture media
- Microbiological media
- Buffer solutions and other biological solutions

Additional Information

No.	6-0007
Manufacturer (Brand)	qpore
EAN	4058072169459
Transport temperature	Room temperature
Color	blue
Volume	250 ml
sterile	No
Material housing	Polystyrene (PS)
Material membrane	Cellulose acetate (CA)
Pore size	0.45 µm
Type connection	Hose olive
Type filter	Bottle top filter
for DM hose inside 1	8-12 mm
for thread/DM	GL 45

