



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## re:Lab qpore® Bottle Top Filter made of CA, sterile, 0.45 µm, 500 ml | 6-0011

Special Price

**€99.00** was

**€199.00**

plus VAT &  
Shipping

### Product Images



## Description

---

**qpore®** offers a comprehensive range of bottle-top filters made of different materials, which are particularly suitable for the filtration of larger sample volumes. All membranes are manufactured under the highest quality standards from the best raw materials.

These **sterile qpore® bottle-top filters with cellulose acetate membrane (CA)** have fast flow rates and very low protein binding, making them ideally suited for filtration of cell culture media. With a GL 45 screw cap and tubing connection, these bottle-top filters are suitable for negative pressure filtration applications directly into the sterile bottle.

### Features:

- Blue color coding
- Quick connector for tubing with inner Ø 8-12 mm
- Graduated filter top
- Sterile, pyrogen free, detergent free
- Bottle-top filters are individually sterile packaged in 24 pieces
- Packaging with cut open aid; thus can be opened without further aids

### Fields of application:

- Sterile filtration of cell culture media
- Microbiological media
- Buffer solutions and other biological solutions

## Additional Information

---

No.	6-0011_RELAB
Manufacturer (Brand)	re:Lab
EAN	4066292691063
Transport temperature	Room temperature
Color	blue
Volume	500 ml
sterile	No
Material membrane	Cellulose acetate (CA)
Pore size	0.45 µm
Type filter	Bottle top filter
for DM hose inside 1	8-12 mm

