



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**qpore® bottle top filter made of PES, sterile, 0.10  
µm, 250 ml, 24 pcs/pack**

**€185.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

**qpore®** offers a comprehensive range of bottle-top filters made of different materials, which are particularly suitable for the filtration of larger sample volumes. All membranes are manufactured under the highest quality standards from the best raw materials.

These **sterile qpore® bottle-top filters with polyethersulfone membrane (PES)** have the fastest flow rates and very low protein binding, making them ideal for filtration of cell culture media. With a GL 45 screw cap and tubing connection, these bottle-top filters are suitable for negative pressure filtration applications directly into the sterile bottle.

### Features:

- Green color coding
- Quick connector for tubing with inner Ø 8-12 mm
- Graduated filter top
- Sterile, pyrogen free, detergent free
- Bottle top filters are individually sterile packed in 24 pieces
- Packaging with cut open aid; thus can be opened without further aids

### Fields of application:

- Sterile filtration of cell culture media
- Microbiological media
- Buffer solutions and other biological solutions

## Additional Information

|                       |                        |
|-----------------------|------------------------|
| No.                   | 6-0037                 |
| Manufacturer (Brand)  | qpore                  |
| EAN                   | 4058072170295          |
| Transport temperature | Room temperature       |
| Color                 | green                  |
| Volume                | 250 ml                 |
| sterile               | Yes                    |
| Material housing      | Polystyrene (PS)       |
| Material membrane     | Polyethersulfone (PES) |
| Pore size             | 0.10 µm                |
| for DM hose inside l  | 8-12 mm                |
| Type connection       | Hose olive             |
| Type filter           | Bottle top filter      |
| for thread/DM         | GL 45                  |

