



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**qpore® bottle-top filter, PVDF, 500 ml, 0.45 µm,  
Ø 90 mm 24 pcs./pack**

**€209.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

**qpore®** offers a comprehensive range of bottle-top filters made of different materials, which are particularly suitable for the filtration of larger sample volumes. All membranes are manufactured under the highest quality standards from the best raw materials.

These **sterile qpore® bottle-top filters with hydrophilic polyvinylidene fluoride (PVDF) membrane** have fastest flow rates. With a GL 45 screw cap and tubing connection, these bottle-top filters are suitable for negative pressure filtration applications directly into the sterile bottle.

### Features:

- Blue color coding
- Quick connector for tubing with inner Ø 8-12 mm
- Graduated filter top
- Sterile, pyrogen free, detergent free
- Bottle-top filters are individually sterile packaged in 24 pieces
- Packaging with cut open aid; thus can be opened without further aids

### Fields of application:

- Sterile filtration of cell culture media
- Microbiological media
- Buffer solutions and other biological solutions

# Additional Information

No.	6-0042
Manufacturer (Brand)	qpore
EAN	4058072168902
Transport temperature	Room temperature
Color	yellow
Volume	500 ml
sterile	Yes
Material housing	Polystyrene (PS)
Material membrane	Polyvinylidene fluoride (PVDF)
Pore size	0.45 µm
for DM hose inside l	8-12 mm
Type connection	Hose olive
Type filter	Bottle top filter
for thread/DM	GL 45

