



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Macherey & Nagel Fafi MN 616 1/4, 18.5 cm

**€41.40**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Pleated filter MN 616 1/4 Diameter: 18.5 cm Pack of 100 pcs.

## Additional Information

---

|                       |                        |
|-----------------------|------------------------|
| No.                   | ML-0940                |
| Dual use indicator    | false                  |
| Manufacturer (Brand)  | Macherey & Nagel       |
| VGKL number           | 104695328              |
| EAN                   | 4058072320973          |
| Manufacturer EAN      | 4046681006142          |
| Gross weight          | 0.2757                 |
| Transport temperature | Room temperature       |
| Net weight            | 0.27570 kg             |
| diameter              | 185 mm                 |
| Country of origin     | Germany                |
| sterile               | No                     |
| Region of origin      | North Rhine-Westphalia |
| Packing width         | 0.19600 m              |
| Packing height        | 0.10600 m              |
| Packing depth         | 0.05000 m              |
| Packaging volume      | 0.00103900 ccm         |
| Customs tariff number | 48232000               |

