



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Ahlstrom pleated filter qualitative 3hw, 185 mm Ø, 100 pcs/pack

**€36.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Pleated filter for testing for contaminants, plant-available nutrients and trace elements in soil samples.

### Specification:

#### Grade 3 hw

Basis weight:65 g/m<sup>2</sup>

Filtration rate:20 s/10 ml

Herzberg:140 s/100 ml

General filtration

## Additional Information

---

No.	E-1417
Manufacturer (Brand)	Ahlstrom
EAN	4058072064617
Transport temperature	Room temperature
diameter	185 mm
Filtration speed	20 sec/10 ml
Basis weight	65 g/m <sup>2</sup>
Heart mountain s/100 ml	220
Surface/ Structure	smooth
Variety Hahnemühle	2555
Variety Munktell	3 hw 15 1002
Variety Whatman	student grade
Type filter	Pleated filter
for soil type	general filtration
for type analysis	qualitative

