



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hahnemühle pleated filter grade 0858, 185 mm, pack of 100 pieces

€42.70
plus VAT &
Shipping

Product Images



Description

Pleated filters (Qualitative) for general laboratory work are recommended for substance identification, clarification of liquids and sample preparation. The filters are made of highly refined cellulose. Granular surface, medium fast filtering.

Applications:

Sample preparations

Clarification of: Alcohols, essences, vinegar, essential oils, extracts/ gelatin, glycerin, paints, lacquers, hair tonics, perfumes, tinctures/ beer wort, spirits, syrup

Use as a protective sheet in filter presses

Additional Information

No.	HN-0073
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4652
EAN	4058072261092
Transport temperature	Room temperature
diameter	185 mm
Basis weight	75 g/m ²
Heart mountain s/100 ml	110
Surface/ Structure	grained
Variety Hahnemühle	858
Variety Munktell	6 49 53
Variety Whatman	183
Type filter	Pleated filter
for soil type	general filtration
for type analysis	qualitative

