



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hahнемühle pleated filter grade 0858, 240 mm, pack of 100 pieces

**€62.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Pleated filters (Qualitative) for general laboratory work are recommended for substance identification, clarification of liquids and sample preparation. The filters are made of highly refined cellulose.  
Granular surface, medium fast filtering.

Applications:

Sample preparations

Clarification of: Alcohols, essences, vinegar, essential oils, extracts/ gelatin, glycerin, paints, lacquers, hair tonics, perfumes, tinctures/  
beer wort, spirits, syrup

Use as a protective sheet in filter presses

## Additional Information

---

No.	HN-0074
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4653
EAN	4058072261108
Transport temperature	Room temperature
diameter	240 mm
Basis weight	75 g/m <sup>2</sup>
for type analysis	qualitative
Heart mountain s/100 ml	110
Surface/ Structure	grained
Variety Hahnemühle	858
Variety Munktell	6 49
Variety Whatman	183
Type filter	Pleated filter
for soil type	general filtration

