



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hahнемühle pleated filter grade 400, 150 mm, pack of 100 pieces

€27.30
plus VAT &
Shipping

Product Images



Description

Pleated filters (Qualitative) for general laboratory work are recommended for substance identification, clarification of liquids and sample preparation. The filters are made of highly refined cellulose.
Smooth surface, medium fast filtering.

Applications:

Sample preparations

Clarification of: Alcohols, essences, vinegar, essential oils, extracts/ gelatin, glycerin, paints, varnishes, hair tonics, perfumes, tinctures/ beer wort, spirits, syrup

Use as a protective sheet in filter presses

Additional Information

No.	HN-0099
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4616
EAN	4058072261191
Transport temperature	Room temperature
diameter	150 mm
Basis weight	61 g/m ²
Heart mountain s/100 ml	200
Surface/ Structure	smooth
Variety Hahnemühle	400
Variety Munktell	4b
Variety Whatman	93
Type filter	Pleated filter
for soil type	general filtration
for type analysis	qualitative

