



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hahнемühle pleated filter grade 400, 185 mm, pack of 100 pieces

**€29.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Pleated filters (Qualitative) for general laboratory work are recommended for substance identification, clarification of liquids and sample preparation. The filters are made of highly refined cellulose.  
Smooth surface, medium fast filtering.

Applications:

Sample preparations

Clarification of: Alcohols, essences, vinegar, essential oils, extracts/ gelatin, glycerin, paints, varnishes, hair tonics, perfumes, tinctures/ beer wort, spirits, syrup

Use as a protective sheet in filter presses

## Additional Information

---

No.	HN-0100
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4617
EAN	4058072261207
Transport temperature	Room temperature
diameter	185 mm
Basis weight	61 g/m <sup>2</sup>
for type analysis	qualitative
Heart mountain s/100 ml	200
Surface/ Structure	smooth
Variety Hahnemühle	400
Variety Munktell	4b
Variety Whatman	93
Type filter	Pleated filter
for soil type	general filtration

