



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hahnemühle pleated filter grade 2555, 320 mm, pack of 100 pieces

€52.40
plus VAT &
Shipping

Product Images



Description

Pleated filters for beer and malt analysis are recommended for brewing analysis methods for filtration and determinations based on MEBAK (Central European Brewing Analysis Commission).

Medium-fast filtration

Ideal for sample preparation and clarification

Suitable for removal of CO₂ and turbidity

Additional Information

No.	HN-0087
Manufacturer (Brand)	Hahnemühle
old neoLab article no.	2-4656
EAN	4058072261139
Transport temperature	Room temperature
diameter	320 mm
Basis weight	75 g/m ²
for type analysis	qualitative
Heart mountain s/100 ml	110
Surface/ Structure	grained
Variety Hahnemühle	2555
Variety Munktell	53
Variety Whatman	183
Type filter	Pleated filter

